Anti-CCP

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Adapted from: Vossenar ER Clin Applied Imm Rev 2004;4:239-262

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Anti-CCP 1 assays are no longer available – when reading papers, be careful to notice what version of anti-CCP assay was used

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Adapted from: Vossenaar ER Clin Applied Imm Rev 2004;4:239-262

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
Adapted from: Vossenar ER and van Venrooij WJ. Clin Applied Imm Rev 2004;4:239-262

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
83 early RA patients

All were former blood donors, so serum samples were available from bloodbank (up to 10 years before first symptoms.)

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ACR Criteria	Odds Ratio	P-value
Morning stiffness >1 hour	2.1	0.108
Arthritis of three or more joints	5.0	0.001
Arthritis of wrist, MCP or PIP	1.2	0.762
Symmetric involvement of joints	6.1	0.002
Rheumatoid nodules	0.003	0.795
Radiographic erosions	8.7	0.001
Positive IgM rheumatoid factor	1.7	0.406
Positive anti-CCP2 antibody	38.6	<0.001

(1 year follow-up)

Presence of anti-CCP in Undifferentiated Arthritis patients accurately predicts future progression to RA.

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American College of Rheumatology (ACR) lineamientos

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Adapted from: Arthritis & Rheumatism 2002;46:328-346

Tratamiento de RA

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Drogas antiinflamatorias, no esteroideas.- Aspirina, Ibuprofeno, Naproxen, etc.

COX-2 inhibitors – e.g., Celecoxib (Celebrex ®)

Esteroides – ej., prednisona, dexametasona.

Disease modifying anti-rheumatic drugs (DMARDs)

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From (09/21/2006): www.glucosamine-arthritis.org/arthritis/RA-Hands.html

Summary (continued)

From: www.rheum.org.uk/health/rheumatoid.html (09/21/2006)

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Source ; AxSYM® Anti-CCP Package insert, immunoassay learning guide


AxSYM® Anti-CCP especificaciones (3)

Comparisión de Métodos	AxSYM Anti-CCP vs. DIASTAT™(Anti-CCP2 assay by Axis-Shield)							
			DIASTAT Anti-CCP		Disease Group		n	% concordance
			POS		NEG			
	AxSYM Anti-CCP	POS	405	1	Non-RA		402	99.5%
		NEG	24	501	RA		529	95.6%
					All samples		931	97.3%
Clinical sensitivity	With the cut-off of 5 U/mL, overall sensitivity was 72.6 %.							
	Disease Group		n		Anti-CCP Positive (n)		% Positive	
	Established RA		331		255		77.0%	
	Early RA		72		48		66.7%	
	Date of Diagnosis unknown		126		81		64.3%	
	All RA Samples		529		384		72.6%	
Clinical specificity	A total 402 specimens, overall specificity in asymptomatic and the various disease groups was 94.5 %.							
					Confirmed Diagnosis			
					POS		NEG	
	AxSYM Anti-CCP	POS	384	22				
		NEG	145	380				

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
**These disease groups gave 100% concordance with the comparison assay (DIASTAT™).*

** Representative data, results in individual laboratories may vary from these data.*

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Disease Group	n	% concordance
Non-RA	402	99.5%
RA	529	95.6%
All samples	931	97.3%

Curve	Area	SE	p	95% CI of Area
DIASTAT	0.870	0.0124	<0.0001	0.846 to 0.895
AxSYM	0.875	0.0119	<0.0001	0.852 to 0.898

Contrast	Difference	p
DIASTAT v AxSYM	-0.005	0.6938

Overall concordance vs. DIASTAT was 97.3 %

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AxSYM[®] Anti-CCP Anti Cyclic Citruillinated Peptide

Ordering Information

AxSYM® Anti-CCP – Ordering Information

Producto		L/N
AxSYM® Anti-CCP Reagent Kit	100 tests	3L91-20
AxSYM® Anti-CCP Calibrators	6 x 4.3 ml	3L91-01
AxSYM® Anti-CCP Controls	2 x 7.0 ml	3L91-10
AxSYM® Anti-CCP Assay Disk	1 disk	3L92-01

AxSYM[®] Anti-CCP Anti Cyclic Citrullinated Peptide

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Anti-mutated citrullinated vimentin (Anti-MCV)

- Orgentec (Germany)
- Very limited data / publications


Anti-Citrullinated Fibrinogen (ACF)

- Biomerieux (VIDAS) - Under development
- (Rumor says development stopped)

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Tabla Comparativa

	<u>AxSYM</u>	<u>Eurodiagnostika Anti-CCP Mk2</u>	<u>ImmunoCAP ELiA</u>
Antigen	CCP2	CCP2	CCP2
Assay technology	Automated MEIA	Manual ELISA microplate	Automated ELiA
Sample Type	Serum/Plasma	Serum/Plasma	Serum/Plasma
Sample Dilution	Automated 1:10	Manual	Automated 1:100
Calibration Stability	Min 14 days	Calibrate every run	1 month
Controls	2 (+/-)	3	2 (+/-)
Sample Volume	150 µl , 83 µl (STAT)	10 µl, manually diluted ad 500 µl	(?)
Primary tube compatible	Yes	No	Yes
Std Curve	6 calibrators	6 calibrators	6 calibrators
Qualitative Option	No	Yes	No
Pack Size	100 tests	96 tests	48 tests
Conjugate	IgG	IgG	IgG
Time to 1 st result	19 mins	2-3 h	2 h
Throughput	Up to 55 tests/hr	96 tests (2-3h)	30 tests/hr (Batch only)
Cutoff value	5 U/ml	25 U/ml	10 U/ml
Clinical specificity	94.5%	97.6%	96.9%
Clinical Sensitivity	72.6%	70%	87.8%
Precision	7.7 - 10.0%CV	6.0 – 12.8%CV	6.6 – 7.0 %CV

Competitive Comparison Table

	<u>AxSYM</u>	<u>Inova QuantaLite CCP</u>	<u>Inova QuantaLite CCP3</u>
Antigen	CCP2	CCP2	"CCP3"
Assay technology	Automated MEIA	Manual ELISA microplate	Manual ELISA microplate
Sample Type	Serum/Plasma	Serum	Serum/Plasma
Sample Dilution	Automated 1:10	Manual	Manual
Calibration Stability	Min 14 days	Calibrate every run	Calibrate every run
Controls	2 (+/-)	2 (+/-)	2 (+/-)
Sample Volume	150 µl , 83 µl (STAT)	5 µl, manually diluted ad 505 µl	5 µl, manually diluted ad 505 µl
Primary tube compatible	Yes	No	No
Std Curve	6 calibrators	5 calibrators	5 calibrators
Qualitative Option	No	No	No
Pack Size	100 tests	96 tests	96 tests
Conjugate	IgG	IgG	IgG
Time to 1 st result	19 mins	90-120 min	90-120 min
Throughput	Up to 55 tests/hr	96 tests (90-120 min)	96 tests (90-120 min)
Cutoff value	5 U/ml	20 U/ml	20 U/ml
Clinical specificity	94.5%	94.8% (98%*)	96%
Clinical Sensitivity	72.6%	76%	74%
Precision	7.7 - 10.0%CV	5 – 11%CV	3 – 8%CV

Competitive Comparison Table

	<u>AxSYM</u>	<u>Euroimmun Anti-CCP</u>	<u>DIASTAT Anti-CCP</u>
Antigen	CCP2	CCP2	CCP2
Assay technology	Automated MEIA	Manual ELISA microplate	Manual ELISA microplate
Sample Type	Serum/Plasma	Serum/Plasma	Serum/Plasma
Sample Dilution	Automated 1:10	Manual	Manual
Calibration Stability	Min 14 days	Calibrate every run	Calibrate every run
Controls	2 (+/-)	?	2 (+/-)
Sample Volume	150 µl , 83 µl (STAT)	?	1-2 µl, use 100 µl 1:100 dilution
Primary tube compatible	Yes	No	No
Std Curve	6 calibrators	5 calibrators	5 calibrators
Qualitative Option	No	?	Yes
Pack Size	100 tests	96 tests	96 tests
Conjugate	IgG	IgG	IgG
Time to 1 st result	19 mins	2 hours +	2 hours +
Throughput	Up to 55 tests/hr	96 tests (over 2 hours)	96 tests (over 2 hours)
Cutoff value	5 U/ml	5 RU/mL	5 U/mL
Clinical specificity	94.5%	98.2%	97%
Clinical Sensitivity	72.6%	78.5%	78%
Precision	7.7 - 10.0%CV	6.3 – 7.2 %	7.6 – 13.6%

Competitive Comparison Table

	<u>AxSYM</u>	<u>Orgentec Anti-MCV</u>
Antigen	CCP2	MCV (mutated citrullinated vimentin)
Assay technology	Automated MEIA	Manual ELISA microplate
Sample Type	Serum/Plasma	Serum
Sample Dilution	Automated 1:10	Manual
Calibration Stability	Min 14 days	Calibrate every run
Controls	2 (+/-)	2 (+/-)
Sample Volume	150 µl , 83 µl (STAT)	10 µl, manually diluted ad 1000 µl
Primary tube compatible	Yes	No
Std Curve	6 calibrators	6 calibrators
Qualitative Option	No	
Pack Size	100 tests	96 tests
Conjugate	IgG	IgG
Time to 1 st result	19 mins	60-90 minutes
Throughput	Up to 55 tests/hr	96 (60-90 min)
Cutoff value	5 U/ml	20 U/ml
Clinical specificity	94.5%	98.7%* *Dejaco et al
Clinical Sensitivity	72.6%	53.7%* *Dejaco et al
Precision	7.7 - 10.0%CV	1.9 – 2.7%CV

AxSYM Anti-CCP Competitors Summary

	<u>ELISA</u>	<u>ImmunoCAP (ELiA)</u>	<u>AxSYM</u>
Antigen	CCP2/CCP3?	CCP2	CCP2
Sample Type	Suero/plasma	Suero/plasma	Suero/plasma
Sample Dilution	Manual	Automated	Automated
Calibration Curve	Cada corrida	Cada 4 semanas	cada 2 a 4 semanas
Software	N/A	Customised	Customised
Sample Dilution	1:100	1:100 (automated)	1:10 (automated)
Sample Volume	10µl	??	150 µl
Primary tube compatible	NO	Yes	Yes
Std Curve	6 Standards	6 Standards	6 calibrators
Qualitative Option	Yes	No	N/A
Pack Size	96 tests	48 tests	100 tests
Conjugate	IgG	IgG	IgG
Time to 1 st result	100-130 mins	120 mins	19 mins
Throughput	30-40 tests/hr	30 tests/hr (Batch only)	53 tests/hr (random)
Bar Coding	No	Yes	Yes
Systems	N/A	~4000 worldwide	>8000 Europe only >16000 worldwide

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